

---

## Product Datasheet

### Anti-Hsp40 (Ydj1) Antibody (orb1149077)

**Description**

Anti-Hsp40 (Ydj1) Mouse Monoclonal

<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Yeast
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IP, WB
<b>Immunogen</b>	Host Species: Mouse Full-length yeast Hsp40 (Ydj1) Antibody Class: IgG1 Preservatives: None.
<b>Target</b>	Hsp40 (Ydj1)
<b>Preservatives</b>	PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P25491</a>
<b>NCBI</b>	<a href="#">NP_014335.1</a>
<b>Entrez</b>	<a href="#">855661</a>
<b>Dilution Range</b>	Immunoblotting: use at 0.5-1ug/ml. A band of ~40 kDa is detected. User should determine optimal concentrations for their application.
<b>Expiration Date</b>	12 months from date of receipt.