
Product Datasheet

Anti-Hsp40 (Ydj1) Antibody (orb1149078)

Description

Anti-HSp40 (Ydj1) Mouse Monoclonal

Species/Host	Mouse
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Host Species: Mouse Full-length yeast Hsp40 (Ydj1) Antibody Class: IgG1 Preservatives: None.
Target	Hsp40 (Ydj1)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P25491
NCBI	NP_014335.1
Entrez	855661
Dilution Range	Immunoblotting: use at 0.5-1ug/ml. A band of ~40 kDa is detected. User should determine optimal concentrations for their application.
Expiration Date	12 months from date of receipt.