

Product Datasheet

Anti-PHAP (CT) Antibody (orb1149101)

Description

Anti-PHAP (CT) Rabbit Polyclonal Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Synthetic peptide corresponding to amino acids near the carboxy terminus of human PHAP (accession no. P39687). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide

Target

PHAP (CT)

Preservatives

PBS, pH 7.4.

Form/Appearance

Liquid

Storage

This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Purity

Purified

Uniprot ID
P39687
NCBI
NP_006296.1
Entrez
8125
Dilution Range

Western Blot: Use at 1ug/ml. Positive control: Raji cell lysate
Immunohistochemistry: use at 10ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.

Expiration Date

12 months from date of receipt.