

Product Datasheet

Anti-PHAP I (CT) Antibody (orb1149102)

Description

Anti-PHAP I (CT) Rabbit Polyclonal Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide corresponding to amino acids at the carboxy terminus of human PHAP (accession no. P39687). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	PHAP I (CT)
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P39687
NCBI	NP_006296.1
Entrez	8125
Dilution Range	Immunoblotting: Use at 1ug/ml. Positive control: Raji cell lysate. Immunohistochemistry: use at 2ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.