
Product Datasheet

Anti-PHAP III Antibody (orb1149104)

Evolution Bioprocessments.
Description

Anti-PHAP III Rabbit Polyclonal Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide corresponding to amino acids near the carboxy terminus of human PHAP III (accession no. NP_112182) Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	PHAP III
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid repeated freezing and thawing.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q9BTT0
NCBI	NP_112182
Entrez	81611
Dilution Range	Western Blot: Use at 1ug/ml. Positive control: A549 cell lysate. Immunohistochemistry: use at 2ug/ml. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.