

Product Datasheet

Anti-ILP-2 Antibody (orb1149125)

Description	Anti-ILP-2 Rabbit Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to aa 2-13 of human ILP-2 (accession no. NP_203127). Host Species: Rabbit Stabilizers: None Preservatives: 0.02% sodium azide
Target	ILP-2
Preservatives	PBS, pH 7.4.
Form/Appearance	Liquid
Storage	This antibody is stable for at least one (1) year at -20C. Avoid multiple freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q96P09
NCBI	NP_203127
Entrez	112401
Dilution Range	Immunoblotting: use at 1ug/ml. Immunocytochemistry: use at 10ug/ml. Positive control: HepG2 or MOLT4 cell lysates. These are recommended concentrations. Enduser should determine optimal concentrations for their applications.
Expiration Date	12 months from date of receipt.