

Product Datasheet

Human Siglec-10 Protein (orb1149297)

Description	Human Siglec-10 Protein, Fc (L117A, L118A) Tag is expressed from human 293 cells (HEK293). It contain...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Siglec-10
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 84.6 kDa. The protein migrates as 100-115 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_149121.2
Purity	95%
Source	Human Siglec-10, Fc Tag (orb1149297) is expressed from human 293 cells (HEK293). It contains AA Met 17 - Thr 546 (Accession # Q96LC7-1).
MW	84.6 kDa
Uniprot ID	Q96LC7
NCBI	NP_149121.2
Expiration Date	6 months from date of receipt.