

Product Datasheet

Human EGFR2 ELISA Kit (orb1149417)

Description	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit
Reactivity	Human
Range	0.5ng/ml-16ng/ml
Tested Applications	ELISA
Target	EGFR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.