

Product Datasheet

Human ADP ELISA Kit (orb1149421)

Description	Human ADP/Acrp30 (Adiponectin) ELISA Kit
Reactivity	Human
Range	0.9375ug/mL-30ug/mL
Tested Applications	ELISA
Target	ADP/Acrp30
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ug/mL
Expiration Date	Please enquire.