

Product Datasheet

Human VCAM-1 ELISA Kit (orb1149423)

Description	Human VCAM-1/CD106 (Vascular Cell Adhesion Molecule 1) ELISA Kit
Reactivity	Human
Range	15ng/ml-480ng/ml
Tested Applications	ELISA
Target	VCAM-1/CD106
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ng/ml
Expiration Date	Please enquire.