

Product Datasheet

Human FNDC5 ELISA Kit (orb1149458)

Description	Human FNDC5 (Fibronectin type III domain-containing protein 5) ELISA Kit
Reactivity	Human
Range	12.5ng/ml-400ng/ml
Tested Applications	ELISA
Target	FNDC5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ng/ml
Expiration Date	Please enquire.