

Product Datasheet

Human VEGFR3 ELISA Kit (orb1149469)

Description	Human VEGFR3/Flt4 (Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit
Reactivity	Human
Range	0.25ng/ml-8ng/ml
Tested Applications	ELISA
Target	VEGFR3/Flt4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.