

Product Datasheet

Human CFHR1 ELISA Kit (orb1149482)

Description	Human CFHR1 (Complement Factor H Related Protein 1) ELISA Kit
Reactivity	Human
Range	25ng/ml-800ng/ml
Tested Applications	ELISA
Target	CFHR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ng/ml
Expiration Date	Please enquire.