

Product Datasheet

Human sCD14 ELISA Kit (orb1149489)

Description	Human sCD14 (Soluble Cluster of Differentiation 14) ELISA Kit
Reactivity	Human
Range	1.875ng/ml-60ng/ml
Tested Applications	ELISA
Target	sCD14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.