

Product Datasheet

Human HSP-27 ELISA Kit (orb1149507)

Description	Human HSP-27/HSPB1 (Heat Shock Protein 27) ELISA Kit
Reactivity	Human
Range	3.75ng/ml-120ng/ml
Tested Applications	ELISA
Target	HSP-27/HSPB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.