

## Product Datasheet

### Human LOX-1 ELISA Kit (orb1149523)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human LOX-1 (Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 10pg/ml-320pg/ml  |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | LOX-1   |
| <b>Storage</b>             | 2-8°C   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Sandwich  |
| <b>Sample Types</b>        | Cells, serum, tissue, saliva, lavage fluid, urine, etc.                         |
| <b>Sensitivity</b>         | 1.0pg/ml  |
| <b>Expiration Date</b>     | Please enquire.   |