

## Product Datasheet

### Human PARC ELISA Kit (orb1149537)

<b>Description</b>	Human PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	2.5ng/ml-80ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	PARC/CCL18
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
<b>Sensitivity</b>	0.1ng/ml
<b>Expiration Date</b>	Please enquire.