

Product Datasheet

Human SP-B ELISA Kit (orb1149539)

Description	Human SP-B (Pulmonary surfactant-associated protein B) ELISA Kit
Reactivity	Human
Range	1.5pg/ml-48pg/ml
Tested Applications	ELISA
Target	SP-B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1pg/ml
Expiration Date	Please enquire.