

## Product Datasheet

### Human LDLR ELISA Kit (orb1149544)

<b>Description</b>	Human LDLR (Low Density Lipoprotein Receptor) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.25mmol/l-8mmol/l
<b>Tested Applications</b>	ELISA
<b>Target</b>	LDLR
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
<b>Sensitivity</b>	0.1mmol/l
<b>Expiration Date</b>	Please enquire.