

## Product Datasheet

### Human VEGFR1 ELISA Kit (orb1149551)

<b>Description</b>	Human VEGFR1/FLT1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	150pg/ml-4800pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	VEGFR1/FLT1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
<b>Sensitivity</b>	10pg/ml
<b>Expiration Date</b>	Please enquire.