

Product Datasheet

Human DBP ELISA Kit (orb1149556)

| | |
|----------------------------|---|
| Description | Human DBP (Vitamin D Binding Protein) ELISA Kit |
| Reactivity | Human |
| Range | 0.25ug/ml-8ug/ml |
| Tested Applications | ELISA |
| Target | DBP |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Cells, serum, tissue, saliva, lavage fluid, urine, etc. |
| Sensitivity | 0.1ug/ml |
| Expiration Date | Please enquire. |