

## Product Datasheet

### Human VEGFR-2 ELISA Kit (orb1149557)

<b>Description</b>	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.625ng/ml-20ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	VEGFR-2/KDR
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
<b>Sensitivity</b>	0.1ng/ml
<b>Expiration Date</b>	Please enquire.