

## Product Datasheet

### Human RBP4 ELISA Kit (orb1149563)

<b>Description</b>	Human RBP4 (Retinol Binding Protein 4) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	2ug/ml-64ug/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	RBP4
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
<b>Sensitivity</b>	0.1ug/ml
<b>Expiration Date</b>	Please enquire.

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)