

Product Datasheet

Human RBP4 ELISA Kit (orb1149563)

Description	Human RBP4 (Retinol Binding Protein 4) ELISA Kit
Reactivity	Human
Range	2ug/ml-64ug/ml
Tested Applications	ELISA
Target	RBP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ug/ml
Expiration Date	Please enquire.