

Product Datasheet

Human PDGF-AB ELISA Kit (orb1149565)

Description	Human PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit
Reactivity	Human
Range	31.25pg/ml-1000pg/ml
Tested Applications	ELISA
Target	PDGF-AB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.