

Product Datasheet

Human HMGB-1 ELISA Kit (orb1149568)

Description	Human HMGB-1 (High mobility group protein B1) ELISA Kit
Reactivity	Human
Range	1250pg/ml-40000pg/ml
Tested Applications	ELISA
Target	HMGB-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10pg/ml
Expiration Date	Please enquire.