

Product Datasheet

Human alpha1-AT ELISA Kit (orb1149580)

Description	Human alpha1-AT (Alpha 1-Antitrypsin) ELISA Kit
Reactivity	Human
Range	125ug/ml-4000ug/ml
Tested Applications	ELISA
Target	α 1-AT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10ug/ml
Expiration Date	Please enquire.