

Product Datasheet

Human TLR-2 ELISA Kit (orb1149581)

Description	Human TLR-2 (Toll-like Receptor 2) ELISA Kit
Reactivity	Human
Range	0.75ng/ml-24ng/ml
Tested Applications	ELISA
Target	TLR-2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.