

Product Datasheet

Human Smad1 ELISA Kit (orb1149583)

Description	Human Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit
Reactivity	Human
Range	0.625ng/ml-20ng/ml
Tested Applications	ELISA
Target	Smad1/MADH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.