

Product Datasheet

Human N-MID-OT ELISA Kit (orb1149586)

Description	Human N-MID-OT (N-MID Osteocalcin) ELISA Kit
Reactivity	Human
Range	0.25ng/ml-8ng/ml
Tested Applications	ELISA
Target	N-MID-OT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.