

Product Datasheet

Human CHEM ELISA Kit (orb1149591)

| | |
|----------------------------|---|
| Description | Human CHEM (Chemerin) ELISA Kit |
| Reactivity | Human |
| Range | 62.5pg/ml-2000pg/ml |
| Tested Applications | ELISA |
| Target | CHEM |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Cells, serum, tissue, saliva, lavage fluid, urine, etc. |
| Sensitivity | 10pg/ml |
| Expiration Date | Please enquire. |