

Product Datasheet

Human AQP-2 ELISA Kit (orb1149613)

| | |
|----------------------------|---|
| Description | Human AQP-2 (Aquaporin 2) ELISA Kit |
| Reactivity | Human |
| Range | 15pg/ml-480pg/ml |
| Tested Applications | ELISA |
| Target | AQP-2 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Cells, serum, tissue, saliva, lavage fluid, urine, etc. |
| Sensitivity | 1.0pg/ml |
| Expiration Date | Please enquire. |