

Product Datasheet

Human AQP-1 ELISA Kit (orb1149614)

Description	Human AQP-1 (Aquaporin 1) ELISA Kit
Reactivity	Human
Range	0.125ng/ml-4ng/ml
Tested Applications	ELISA
Target	AQP-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.