

Product Datasheet

Human ApoC2 ELISA Kit (orb1149619)

Description	Human ApoC2 (Apolipoprotein C2) ELISA Kit
Reactivity	Human
Range	0.375mg/dl-12mg/dl
Tested Applications	ELISA
Target	ApoC2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1mg/dl
Expiration Date	Please enquire.