

Product Datasheet

Human BRAK ELISA Kit (orb1149627)

Description	Human BRAK/CXCL14 (Breast and Kidney Expressed Chemokine) ELISA Kit
Reactivity	Human
Range	31.25pg/ml-1000pg/ml
Tested Applications	ELISA
Target	BRAK/CXCL14
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.