

Product Datasheet

Human IGFBP-4 ELISA Kit (orb1149635)

Description	Human IGFBP-4 (Insulin-like Growth Factor Binding Protein 4) ELISA Kit
Reactivity	Human
Range	5ng/ml-160ng/ml
Tested Applications	ELISA
Target	IGFBP-4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0ng/ml
Expiration Date	Please enquire.