

Product Datasheet

Human IGF1R ELISA Kit (orb1149639)

Description	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit
Reactivity	Human
Range	6.25ng/ml-200ng/ml
Tested Applications	ELISA
Target	IGF1R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0g/ml
Expiration Date	Please enquire.