

Product Datasheet

Human ANGPTL3 ELISA Kit (orb1149644)

Description	Human ANGPTL3 (Angiopietin Like Protein 3) ELISA Kit
Reactivity	Human
Range	2.5ng/ml-80ng/ml
Tested Applications	ELISA
Target	ANGPTL3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.