

Product Datasheet

Human IL-1beta ELISA Kit (orb1149673)

Description	Human IL-1beta (Interleukin 1 Beta) ELISA Kit
Reactivity	Human
Range	2.5pg/ml-80pg/ml
Tested Applications	ELISA
Target	IL-1 β
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1pg/ml
Expiration Date	Please enquire.