

Product Datasheet

Human EGFR ELISA Kit (orb1149676)

Description	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit
Reactivity	Human
Range	0.625ng/ml-200ng/ml
Tested Applications	ELISA
Target	EGFR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.