

## Product Datasheet

### Human MIP-1alpha ELISA Kit (orb1149691)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human MIP-1alpha (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 5ng/ml-160ng/ml  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | MIP-1 $\alpha$   |
| <b>Storage</b>             | 2-8°C  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Sandwich   |
| <b>Sample Types</b>        | Cells, serum, tissue, saliva, lavage fluid, urine, etc.              |
| <b>Sensitivity</b>         | 1.0ng/ml   |
| <b>Expiration Date</b>     | Please enquire.  |