

Product Datasheet

Human BMP-4 ELISA Kit (orb1149693)

Description	Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit
Reactivity	Human
Range	6.25pg/ml-200pg/ml
Tested Applications	ELISA
Target	BMP-4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.