

Product Datasheet

Human BAFF ELISA Kit (orb1149696)

Description	Human BAFF/CD257 (B-Cell Activating Factor) ELISA Kit
Reactivity	Human
Range	37.5pg/ml-1200pg/ml
Tested Applications	ELISA
Target	BAFF/CD257
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.