

Product Datasheet

Human IGFBP-3 ELISA Kit (orb1149721)

Description	Human IGFBP-3 (Insulin-like Growth Factor Binding Protein 3) ELISA Kit
Reactivity	Human
Range	3.75ng/ml-120ng/ml
Tested Applications	ELISA
Target	IGFBP-3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	0.1ng/ml
Expiration Date	Please enquire.