

Product Datasheet

Human GM-CSF ELISA Kit (orb1149725)

Description	Human GM-CSF (Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit
Reactivity	Human
Range	20pg/ml-640pg/ml
Tested Applications	ELISA
Target	GM-CSF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	1.0pg/ml
Expiration Date	Please enquire.