

Product Datasheet

Human GDF15 ELISA Kit (orb1149726)

Description	Human GDF15 (Growth Differentiation Factor 15) ELISA Kit
Reactivity	Human
Range	50pg/ml-1600pg/ml
Tested Applications	ELISA
Target	GDF15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10pg/ml
Expiration Date	Please enquire.