

Product Datasheet

Human Flt3L ELISA Kit (orb1149729)

Description	Human Flt3L (FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit
Reactivity	Human
Range	62.5pg/ml-2000pg/ml
Tested Applications	ELISA
Target	Flt3L
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Cells, serum, tissue, saliva, lavage fluid, urine, etc.
Sensitivity	10pg/ml
Expiration Date	Please enquire.