

Product Datasheet

Symmetric DiMethyl-Histone H3-R8 antibody (orb1149740)

Description	Rabbit polyclonal antibody to Symmetric DiMethyl-Histone H3-R8
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	A synthetic symmetric dimethylated peptide around R8 of human histone H3(NP_003520.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17KDa
Uniprot ID	P68431
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.