

## Product Datasheet

### SHC1 (Phospho-Y427) antibody (orb1149762)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to SHC1 (Phospho-Y427)   |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | A synthetic phosphorylated peptide around Y427 of human SHC1(NP_892113.4).  |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 50KDa/55KDa   |
| <b>Uniprot ID</b>          | <b>P29353</b>   |
| <b>Dilution Range</b>      | WB 1:500 - 1:2000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |