

## Product Datasheet

# Asymmetric DiMethyl-Histone H4-R3 antibody (orb1149794)

<b>Description</b>	Rabbit polyclonal antibody to Asymmetric DiMethyl-Histone H4-R3
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	A synthetic peptide of human Asymmetric DiMethyl-Histone H4-R3.
<b>Preservatives</b>	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>Uniprot ID</b>	<b>P62805</b>
<b>Entrez</b>	<b>8359</b>
<b>Expiration Date</b>	12 months from date of receipt.